

......

## Aegis Remote Monitoring IoT System

## **TO ORDER**

Call 1-800-796-2344 ask for Elizabeth Hess OR send Email: orders@aegisfridge.com For support send email to support@aegisfridge.com

## HOW IT WORKS REMOTE MONITORING SYSTEM

Aegis Remote Monitoring IoT System consists of (i) Sensors, (ii) Communication Network, (iii) Cloud and (iv) Dashboard. The sensors are modular in design with a digital interface for different types of sensors probes such as T, RH, CO2, Differential pressure and more. The sensor platform is modular which makes easy addition of sensors, including any off-the-shelf sensor with customization. Several diagnostic tools have been developed to measure the sensor health remotely allowing the solution to be implemented on a large scale.



- 1. Monitored Refrigerator or Freezer
- 2. IoT Cloud with Ingestion, Data Processing, Storage and Analytics
- 3. End User Dashboard for Monitoring

The sensor communicates using WiFi Protocols (IEEE 802.11b/g/n). A keep feature of our WiFi communication is that it supports a myriad of encryption including WEP, WPA/WPA2, PEAP, EAP-TLS. Therefore, our sensors can easily operate in any corporate or enterprise network. On board data storage with the store-and-forward feature prevents any data loss in case of network problems. Use of WiFi eliminates the need for any gateway device. The sensors can communicate directly to the Cloud. Fig 1. shows the data flow from the sensor to the cloud.





Fig. 1

The IoT cloud platform is unique because it has an **ingestion engine that supports virtually unlimited number of sensors across a distributed network over many locations**. The key elements of the IoT cloud platform is depicted in Fig 2. An additional benefit of this platform is the possibility for users to develop prediction models for a particular use case using the AI and machine learning modules.

DASHBO	DARD	MANAGEN	IENT	ALARM		REPORTS	ADMIN	3.7.7.8	octavche@gmail.co
Sensors	Assets		Hards		<b>A</b> rativo	0.4.4	1	Name: Sensor_00202562 EDIT	2 ×
Sensors	Battery	Sending S Data	Weak Signal 31/195	Low Uptime 173/195	Active Alarms 2/195	Out-of- date firmware 194/195	Has stored data 10/195	ID: 00202562 Product Number: External Pro Probe - XTEMP_2000_0001 Alarm Prome. SR P   Asset: None Location: San Rafael Clinic Lab Alarm enabled: true   Group: San Rafael Clinic Coalition: Marin Community Clinics Repeat alarming: n   Sampling: 300s Upload Period: 300s Automatically close   Code Version: 7.3.39 Code Version: 7.3.39 Automatically close	ever
Search sense	ors							Last Seen: 2m Last Data: 🗥 2m Link Quality -53 dBm Battery	0.00 shareiza
Name		ID	T1	T2	Last Seen	<b>۱</b> 个		Last Seen: 2m Last Data: 2m Link Quality -53 dBm Battery Voltage	3.99 charging V ~90%
UG115F	Room	00200115	22.92	N/A	Now	🛛 🌲 🛡	1	temp1: 5.36 °C 1.90 °C 8.10 °C	
🔥 Mason N	lo Probesz	00200855	N/A	N/A	Now	or 🖄 🗢	۲ <u>۵</u>		
🔥 Sensor_C	0202091	00202091	28.94	N/A	Now	📀 🔌 🗬	۵.	Graph Data Open Alarms Closed Alarms Reports Health	3
🔥 Contacto	orCoilB Temp	00202232	57.50	N/A	Now	O 🌲 🗣	۲ <u>۵</u>		
Sensor_C	0202566	00202566	3.99	N/A	Now	O 🌲 🛡	P 🖬	Type   Date   From   To   Show datapoin     temp1 *C	ts Show battery & signal
Sensor_C	0203050	00203050	28.26	N/A	Now	<b>I</b>	2		Ξ
Bensor_0	0202615	00202615	5.64	N/A	1m	O 🌲 🛡	<b>D</b>	Upper threshold (8.1°C) Safe area	<u>&amp;°C</u>
U Sensor_C	0202559	00202559	-21.92	N/A	1m	🗢 🌲 🛡	۵.	Max	6°C
U Sensor_C	0202562	00202562	5.36	N/A	1m	o 🇯 🗢	â		
() Contacto	orCoilA temp	00202234	57.20	N/A	2m	<b>&gt;</b> 🔌 🛡	D D		4°C
Sensor_C	0202611	00202611	4.10	N/A	2m	<b>S</b>	ù -	Lower threshold (1.9°C)	2°C
Sensor_C	0202605	00202605	-20.16	N/A	2m	<b>S</b>	P Q		0°C
									Fig. 3

- 1. Sensor List
- 2. Selected Sensor info and sensor health
- 3. Sensor reading (Graph, Data, Alarms, Reports, Health)

From the cloud the sensors data is moved to a **dashboard for data visualization**, **alerting and reporting**. Alerting engine includes emails and text messages with escalations. Custom reporting features, for any particular domain, can be created on the dashboard. The dashboard can be accessed from a desktop computer or a phone or a tablet. Fig 3. is an example of the dashboard with the sensor list on the left panel and sensor details on the right panel.

Another unique feature of the dashboard is to track annual sensor certifications and validation procedures. In addition the ability to visualize assets and its maintenance procedures is very useful. Custom workflows can be implemented.

